

# NITRILE

## POWDER-FREE GLOVES



Made of synthetic rubber, contain no latex. Excellent puncture resistant and excellent resistance to wear and tears. High degree of flexibility and superior resistance to solvent / many types of chemicals. Thin nitrile offers excellent strength and great comfort.

### PHYSICAL PROPERTIES

	Force at Break EN455-2	Tensile Strength ASTM D6319	Elongation ASTM D6319	Elongation ISO 11193-1	Force at Break ISO 11193-1
Before Aging	≥6N	≥14Mpa	≥500%	≥500%	≥7N
After Aging	≥6N	≥14Mpa	≥400%	≥400%	≥7N

### DIMENSION & TOLERANCES

Extra Small (mm)	: 75 ±5mm
Small (mm)	: 85 ±5mm
Medium (mm)	: 95 ±5mm
Large (mm)	: 105 ±5mm
Extra Large (mm)	: 115 ±5mm
Cuff (25mm from the bead)	: Min 0.04mm
Palm (Center of the palm)	: Min 0.05mm
Finger (13mm from tip)	: Min 0.07mm
Length (mm)	: Min 240mm

### PRODUCT SPECIFICATIONS

Regulatory Compliance	: FDA 510K, MDD 93/42/EEC, Regulation (EU) 2016/425, ISO 13485, ISO 9001
Standard	: ASTM D6319, EN 455, ISO 11193-1, EN374, ISO 16604
Classification	: Class I (FDA), Class I (MDD 93/42/EEC), Class III (PPE)
Material	: Nitrile Butadiene Rubber
Color	: Blue/ White/ Other Coloured Gloves
Protein Content	: N/A
Powder Residue Content	: Less than 2 mg/glove

Storage condition (recommended): Cool & dry place with temperature below 40°C.



博利手套  
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia. Fax: (605) 891 1417

www.professionallatex.com

# 丁腈无粉 乳胶检查手套



合成橡胶制成, 不含天然乳胶蛋白质。比其他手套更优异的穿刺抗性和更强的耐磨性。它具有卓越的触感和灵活性 / 高度耐抗多种化学品。轻薄的丁腈手套有高抗拉力强度, 为使用者提供了优质柔软与舒适感。

## 物理性质

	断点强力 EN455-2	抗拉力强度 ASTM D6319	伸长率 ASTM D6319	伸长率 ISO 11193-1	断点强力 ISO 11193-1
老化前	≥6N	≥14Mpa	≥500%	≥500%	≥7N
老化后	≥6N	≥14Mpa	≥400%	≥400%	≥7N


## 尺寸和公差

特小 (毫米)	: 75 ± 5毫米
小 (毫米)	: 85 ± 5毫米
中 (毫米)	: 95 ± 5毫米
大 (毫米)	: 105 ± 5毫米
特大 (毫米)	: 115 ± 5毫米
袖口 (25毫米从袖口)	: 不低于 0.04毫米
手掌部 (手掌中心)	: 不低于 0.05毫米
手指部 (13毫米从尖)	: 不低于 0.07毫米
长度 (毫米)	: 不低于 240毫米

## 产品规格表

质量控管	: FDA 510K, MDD 93/42/EEC, Regulation (EU) 2016/425, ISO 13485, ISO 9001
产品质量标准	: ASTM D6319, EN 455, ISO 11193-1, EN374, ISO 16604
分类	: Class I (FDA), Class I (MDD 93/42/EEC), Class III (PPE)
材料	: 丁腈橡胶
颜色	: 蓝/白/其他颜色手套 (根据客户要求)
蛋白质含量	: N/A
粉末残留量	: 少于 2 mg/glove

存储条件 (推荐): 凉爽干燥的地方, 温度低于 40° C.



博利手套  
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. [1038089-A]

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia. Fax: (605) 891 1417

www.professionallatex.com