

# CHLORINATED EXAMINATION GLOVES

## NATURAL RUBBER GLOVE

The natural rubber glove provides superior barrier protection with excellent tensile strength and tear resistant. This gloves are manufactured under clean and safe environment with lowest allowable protein claim for natural latex glove.



### PHYSICAL PROPERTIES

	Before Aging	After Aging
<b>Force at Break</b> EN455-2	≥6N	≥6N
<b>Tensile Strength</b> ASTM D3578	≥18MPa	≥14MPa
<b>Elongation</b> ASTM D3578	≥650%	≥500%
<b>Elongation</b>	≥650%	≥500%
<b>Force at Break</b>	≥7N	≥6N

### DIMENSION & TOLERANCES

Extra Small (mm)	: 75 ±5mm	Cuff (25mm from the bead)	: Min 0.08mm
Small (mm)	: 85 ±5mm	Palm (Center of the palm)	: Min 0.10mm
Medium (mm)	: 95 ±5mm	Finger (13mm from tip)	: Min 0.12mm
Large (mm)	: 105 ±5mm	Length (mm)	: Min 240mm
Extra Large (mm)	: 115 ±5mm		

### PRODUCT SPECIFICATIONS

Regulatory Compliance	: FDA 510K, MDR 2017/745, Regulation (EU) 2016/425, ISO 13485, ISO 9001
Standard	: ASTM D3578, EN 455, EN 14362-1:2017, EN374, ISO 16604
Classification	: Class I (FDA), Class I (MDR 2017/745), Cat III (PPE)
Material	: Natural Rubber Latex
Color	: Yellow (DC) or Light Yellow (AC)
Protein Content	: Less than 50 µg/dm <sup>2</sup>
Powder Residue Content	: Less than 2 mg/glove

Storage condition (recommended): Cool & dry place with temperature below 40°C.



博利手套  
Professional Latex

+605 891 6650 [marketing@professionallatex.com](mailto:marketing@professionallatex.com)

PROFESSIONAL LATEX SDN. BHD. (1038089-A)  
Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia.  
Fax: (605) 891 1417

[www.professionallatex.com](http://www.professionallatex.com)

# 无粉氯洗无灭菌 乳胶检查手套

## 天然乳胶手套

天然乳胶手套提供优越的保护屏障和优越的抗拉力强度。此手套是在清洁和安全的环境下生产的，低含量蛋白质及杂物。



## 物理性质

	老化前	老化后
<b>断点强力</b> EN455-2	≥6N	≥6N
<b>抗拉力强度</b> ASTM D3578	≥18MPa	≥14MPa
<b>伸长率</b> ASTM D3578	≥650%	≥500%
<b>伸长率</b>	≥650%	≥500%
<b>断点强力</b>	≥7N	≥6N

## 尺寸和公差

特小 (毫米)	: 75 ±5毫米	袖口 (25毫米从袖口)	: 不低于 0.08毫米
小 (毫米)	: 85 ±5毫米	手掌部 (手掌中心)	: 不低于 0.10毫米
中 (毫米)	: 95 ±5毫米	手指部 (13毫米从尖)	: 不低于 0.12毫米
大 (毫米)	: 105 ±5毫米	长度 (毫米)	: 不低于 240毫米
特大 (毫米)	: 115 ±5毫米		

## 产品规格表

质量控管	: FDA 510K, MDR 2017/745, Regulation (EU) 2016/425, ISO 13485, ISO 9001
产品质量标准	: ASTM D3578, EN 455, EN 14362-1:2017, EN374, ISO 16604
分类	: Class I (FDA), Class I (MDR 2017/745), Cat III (PPE)
材料	: 天然胶乳
颜色	: 黄色 (DC) 或浅黄色 (AC)
蛋白质含量	: 少于 50µg/dm <sup>2</sup>
粉末残留量	: 少于 2 mg/glove

存储条件 (推荐): 凉爽干燥的地方, 温度低于 40°C.



博利手套  
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia.

Fax: (605) 891 1417

www.professionallatex.com