

NITRILE POWDER-FREE GLOVES

MADE OF SYNTHETIC RUBBER, LATEX FREE

Excellent puncture resistant and excellent resistance to wear and tears. High degree of flexibility and superior resistance to solvent / many types of chemicals. Thin nitrile offers excellent strength and great comfort.



PHYSICAL PROPERTIES

	Before Aging	After Aging
Force at Break EN455-2	≥6N	≥6N
Tensile Strength ASTM D6319	≥14MPa	≥14MPa
Elongation ASTM D6319	≥500%	≥400%
Elongation	≥500%	≥400%
Force at Break	≥7N	≥7N

DIMENSION & TOLERANCES

Extra Small (mm)	: 75 ±5mm	Cuff (25mm from the bead)	: Min 0.04mm
Small (mm)	: 85 ±5mm	Palm (Center of the palm)	: Min 0.05mm
Medium (mm)	: 95 ±5mm	Finger (13mm from tip)	: Min 0.06mm
Large (mm)	: 105 ±5mm	Length (mm)	: Min 240mm
Extra Large (mm)	: 115 ±5mm		

PRODUCT SPECIFICATIONS

Regulatory Compliance	: FDA 510K, MDR 2017/745, Regulation (EU) 2016/425, ISO 13485, ISO 9001
Standard	: ASTM D6319, EN 455, EN 14362-1:2017, EN374, ISO 16604
Classification	: Class I (FDA), Class I (MDR 2017/745), Cat III (PPE)
Material	: Nitrile Butadiene Rubber
Color	: Blue/ White/ Other Coloured Gloves
Protein Content	: N/A
Powder Residue Content	: Less than 2 mg/glove

Storage condition (recommended): Cool & dry place with temperature below 40°C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)
Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia.
Fax: (605) 891 1417

www.professionallatex.com

丁腈无粉 乳胶检查手套

合成橡胶制成，不含天然乳胶蛋白质

比其他手套更优异的穿刺抗性和更强的耐磨性。它具有卓越的触感和灵活性 / 高度耐抗多种化学品。轻薄的丁腈手套有高抗拉力强度，为使用者提供了优质柔软与舒适感。



物理性质

	老化前	老化后
断点强力 EN455-2	≥6N	≥6N
抗拉力强度 ASTM D6319	≥14MPa	≥14MPa
伸长率 ASTM D6319	≥500%	≥400%
伸长率	≥500%	≥400%
断点强力	≥7N	≥7N

尺寸和公差

特小 (毫米)	: 75 ±5毫米	袖口 (25毫米从袖口)	: 不低于 0.04毫米
小 (毫米)	: 85 ±5毫米	手掌部 (手掌中心)	: 不低于 0.05毫米
中 (毫米)	: 95 ±5毫米	手指部 (13毫米从尖)	: 不低于 0.06毫米
大 (毫米)	: 105 ±5毫米	长度 (毫米)	: 不低于 240毫米
特大 (毫米)	: 115 ±5毫米		

产品规格表

质量控管	: FDA 510K, MDR 2017/745, Regulation (EU) 2016/425, ISO 13485, ISO 9001
产品质量标准	: ASTM D6319, EN 455, EN 14362-1:2017, EN374, ISO 16604
分类	: Class I (FDA), Class I (MDR 2017/745), Cat III (PPE)
材料	: 丁腈橡胶
颜色	: 蓝/白/其他颜色手套 (根据客户要求)
蛋白质含量	: N/A
粉末残留量	: 少于 2 mg/glove

存储条件 (推荐): 凉爽干燥的地方, 温度低于 40°C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)
Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia.
Fax: (605) 891 1417

www.professionallatex.com