

POLYMER COATED EXAMINATION GLOVES

Our Polymer coated latex gloves with pleasant odour, superior aging characteristic and fully textured finish for excellent grip.



PHYSICAL PROPERTIES

	Before Aging	After Aging
Force at Break EN455-2	≥6N	≥6N
Tensile Strength ASTM D3578	≥18MPa	≥14MPa
Elongation ASTM D3578	≥650%	≥500%
Elongation	≥650%	≥500%
Force at Break	≥7N	≥6N

DIMENSION & TOLERANCES

Extra Small (mm)	: 75 ±5mm	Cuff (25mm from the bead)	: Min 0.08mm
Small (mm)	: 85 ±5mm	Palm (Center of the palm)	: Min 0.10mm
Medium (mm)	: 95 ±5mm	Finger (13mm from tip)	: Min 0.12mm
Large (mm)	: 105 ±5mm	Length (mm)	: Min 240mm
Extra Large (mm)	: 115 ±5mm		

PRODUCT SPECIFICATIONS

Regulatory Compliance	: FDA 510K, MDR 2017/745, Regulation (EU) 2016/425, ISO 13485, ISO 9001
Standard	: ASTM D3578, EN 455, EN 14362-1:2017, EN374, ISO 16604
Classification	: Class I (FDA), Class I (MDR 2017/745), Cat III (PPE)
Material	: Natural Rubber Latex
Color	: Natural White
Protein Content	: Less than 50 µg/dm ²
Powder Residue Content	: Less than 2 mg/glove

Storage condition (recommended): Cool & dry place with temperature below 40°C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)
Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia.
Fax: (605) 891 1417

www.professionallatex.com

无粉涂层无灭菌 乳胶检查手套

我们的聚合物层手套无异味，有优良的老化特性，充分的质感处理让手套抓力更强。



无异味



优良的
老化特性



物理性质

	老化前	老化后
断点强力 EN455-2	≥6N	≥6N
抗拉力强度 ASTM D3578	≥18MPa	≥14MPa
伸长率 ASTM D3578	≥650%	≥500%
伸长率	≥650%	≥500%
断点强力	≥7N	≥6N

尺寸和公差

特小 (毫米)	: 75 ±5毫米	袖口 (25毫米从袖口)	: 不低于 0.08毫米
小 (毫米)	: 85 ±5毫米	手掌部 (手掌中心)	: 不低于 0.10毫米
中 (毫米)	: 95 ±5毫米	手指部 (13毫米从尖)	: 不低于 0.12毫米
大 (毫米)	: 105 ±5毫米	长度 (毫米)	: 不低于 240毫米
特大 (毫米)	: 115 ±5毫米		

产品规格表

质量控管	: FDA 510K, MDR 2017/745, Regulation (EU) 2016/425, ISO 13485, ISO 9001
产品质量标准	: ASTM D3578, EN 455, EN 14362-1:2017, EN374, ISO 16604
分类	: Class I (FDA), Class I (MDR 2017/745), Cat III (PPE)
材料	: 天然乳胶
颜色	: 白色
蛋白质含量	: 少于 50µg/dm ²
粉末残留量	: 少于 2 mg/glove

存储条件 (推荐): 凉爽干燥的地方, 温度低于 40°C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia.

Fax: (605) 891 1417

www.professionallatex.com