

CHLORINATED EXAMINATION GLOVES



The natural rubber glove provides superior barrier protection with excellent tensile strength and tear resistant. This gloves are manufactured under clean and safe environment with lowest allowable protein claim for natural latex glove.

PHYSICAL PROPERTIES

	Force at Break EN455-2	Tensile Strength ASTM D3578	Elongation ASTM D3578	Elongation ISO 11193-1	Force at Break ISO 11193-1
Before Aging	≥6N	≥18Mpa	≥650%	≥650%	≥7N
After Aging	≥6N	≥14Mpa	≥500%	≥500%	≥6N


DIMENSION & TOLERANCES

Extra Small (mm)	: 75 ±5mm
Small (mm)	: 85 ±5mm
Medium (mm)	: 95 ±5mm
Large (mm)	: 105 ±5mm
Extra Large (mm)	: 115 ±5mm
Cuff (25mm from the bead)	: Min 0.07mm
Palm (Center of the palm)	: Min 0.09mm
Finger (13mm from tip)	: Min 0.11mm
Length (mm)	: Min 240mm

PRODUCT SPECIFICATIONS

Regulatory Compliance	: FDA 510K, MDD 93/42/EEC, Regulation (EU) 2016/425, ISO 13485, ISO 9001
Standard	: ASTM D3578, EN 455, ISO 11193-1, EN374, ISO 16604
Classification	: Class I (FDA), Class I (MDD 93/42/EEC), Class III (PPE)
Material	: Natural Rubber Latex
Color	: Natural Light Yellow
Protein Content	: Less than 50 µg/dm ²
Powder Residue Content	: Less than 2 mg/glove

Storage condition (recommended): Cool & dry place with temperature below 40°C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia. Fax: (605) 891 1417

www.professionallatex.com

无粉无灭菌 乳胶检查手套



天然乳胶手套提供优越的保护屏障和优越的抗拉力强度。此手套是在清洁和安全的环境下生产的，低含量蛋白质及杂物。

物理性质

	断点强力 EN455-2	抗拉力强度 ASTM D3578	伸长率 ASTM D3578	伸长率 ISO 11193-1	断点强力 ISO 11193-1
--	-----------------	---------------------	-------------------	--------------------	---------------------

老化前	≥6N	≥18Mpa	≥650%	≥650%	≥7N
-----	-----	--------	-------	-------	-----

老化后	≥6N	≥14Mpa	≥500%	≥500%	≥6N
-----	-----	--------	-------	-------	-----

尺寸和公差

特小 (毫米)	: 75 ± 5毫米
小 (毫米)	: 85 ± 5毫米
中 (毫米)	: 95 ± 5毫米
大 (毫米)	: 105 ± 5毫米
特大 (毫米)	: 115 ± 5毫米
袖口 (25毫米从袖口)	: 不低于 0.07毫米
手掌部 (手掌中心)	: 不低于 0.09毫米
手指部 (13毫米从尖)	: 不低于 0.11毫米
长度 (毫米)	: 不低于 240毫米

产品规格表

质量控管	: FDA 510K, MDD 93/42/EEC, Regulation (EU) 2016/425, ISO 13485, ISO 9001
产品质量标准	: ASTM D3578, EN 455, ISO 11193-1, EN374, ISO 16604
分类	: Class I (FDA), Class I (MDD 93/42/EEC), Class III (PPE)
材料	: 天然胶乳
颜色	: 浅黄色
蛋白质含量	: 少于 50 µg/dm ²
粉末残留量	: 少于 2 mg/glove

存储条件 (推荐): 凉爽干燥的地方, 温度低于 40° C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. [1038089-A]

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia. Fax: (605) 891 1417

www.professionallatex.com