

POLYMER COATED EXAMINATION GLOVES



Our Polymer coated latex gloves with pleasant odour, superior aging characteristic and fully textured finish for excellent grip.

PHYSICAL PROPERTIES

	Force at Break EN455-2	Tensile Strength ASTM D3578	Elongation ASTM D3578	Elongation ISO 11193-1	Force at Break ISO 11193-1
Before Aging	≥6N	≥18Mpa	≥650%	≥650%	≥7N
After Aging	≥6N	≥14Mpa	≥500%	≥500%	≥6N

DIMENSION & TOLERANCES

Extra Small (mm)	: 75 ±5mm
Small (mm)	: 85 ±5mm
Medium (mm)	: 95 ±5mm
Large (mm)	: 105 ±5mm
Extra Large (mm)	: 115 ±5mm
Cuff (25mm from the bead)	: Min 0.07mm
Palm (Center of the palm)	: Min 0.11mm
Finger (13mm from tip)	: Min 0.13mm
Length (mm)	: Min 240mm

PRODUCT SPECIFICATIONS

Regulatory Compliance	: FDA 510K, MDD 93/42/EEC, Regulation (EU) 2016/425, ISO 13485, ISO 9001
Standard	: ASTM D3578, EN 455, ISO 11193-1, EN374, ISO 16604
Classification	: Class I (FDA), Class I (MDD 93/42/EEC), Class III (PPE)
Material	: Natural Rubber Latex
Color	: Natural White
Protein Content	: Less than 50 µg/dm ²
Powder Residue Content	: Less than 2 mg/glove

Storage condition (recommended): Cool & dry place with temperature below 40°C.



博利手套
Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia. Fax: (605) 891 1417

www.professionallatex.com

无粉无灭菌 乳胶检查手套



我们的聚合物层手套无异味，有优良的老化特性，充分的质感处理让手套抓力更强。

物理性质

	断点强力 EN455-2	抗拉力强度 ASTM D3578	伸长率 ASTM D3578	伸长率 ISO 11193-1	断点强力 ISO 11193-1
老化前	≥6N	≥18Mpa	≥650%	≥650%	≥7N
老化后	≥6N	≥14Mpa	≥500%	≥500%	≥6N


尺寸和公差

特小 (毫米)	: 75 ± 5毫米
小 (毫米)	: 85 ± 5毫米
中 (毫米)	: 95 ± 5毫米
大 (毫米)	: 105 ± 5毫米
特大 (毫米)	: 115 ± 5毫米
袖口 (25毫米从袖口)	: 不低于 0.07毫米
手掌部 (手掌中心)	: 不低于 0.11毫米
手指部 (13毫米从尖)	: 不低于 0.13毫米
长度 (毫米)	: 不低于 240毫米

产品规格表

质量控管	: FDA 510K, MDD 93/42/EEC, Regulation (EU) 2016/425, ISO 13485, ISO 9001
产品质量标准	: ASTM D3578, EN 455, ISO 11193-1, EN374, ISO 16604
分类	: Class I (FDA), Class I (MDD 93/42/EEC), Class III (PPE)
材料	: 天然乳胶
颜色	: 白色
蛋白质含量	: 少于 50 µg/dm ²
粉末残留量	: 少于 2 mg/glove

存储条件 (推荐): 凉爽干燥的地方, 温度低于 40° C.



博利手套

Professional Latex

+605 891 6650 marketing@professionallatex.com

PROFESSIONAL LATEX SDN. BHD. (1038089-A)

Lot. 52, Jalan Logam 2, Kamunting Raya Industrial, 34600, Kamunting, Perak, Malaysia. Fax: (605) 891 1417

www.professionallatex.com