

## NATURAL RUBBER

# CHLORINATED LATEX GLOVES


These powder-free, biodegradable natural rubber gloves offer enhanced tactile sensitivity, comfort, and excellent barrier protection against chemicals, viruses, and hazardous substances. Compliant with FDA food handling requirements, they come in various textured surfaces and feature an ambidextrous design with a beaded cuff.



Powder Free



Ambidextrous



Chlorinated  
(Beaded: 240mm)



Weight Range  
5.5g ~ 8.0g

### Standards & Certifications



### Available in

Standard Chlorinated Latex Gloves



Chlorinated Latex EXAM Gloves



### Colour

Light Yellow

### Application



General Use



Medical



General Lab Work



Healthcare



Life Sciences



Dental




Chemical



Semiconductor



Scan the QR Code to  
Download Our Latest Catalogue



# LATEX GLOVES

Polymer Latex Gloves	Chlorinated Latex Gloves	Chlorinated Latex EXAM Gloves	Double Chlorinated Gloves
<ul style="list-style-type: none"> <li>Complies to FDA food handling requirements</li> </ul>	<ul style="list-style-type: none"> <li>Complies to FDA food handling requirements</li> <li>Good barrier protection against chemicals, virus and dangerous substances</li> </ul>	<ul style="list-style-type: none"> <li>Good barrier protection against chemicals, virus and dangerous substances</li> <li>Complies to EU MDR requirements</li> </ul>	-

- Powder Free
- Latex gloves designed to provide enhanced tactile sensitivity
- Soft, elastic latex material provide excellent comfort
- Available in finger textured, smooth, and textured surface
- Ambidextrous, beaded cuff



## APPLICATIONS

<ul style="list-style-type: none"> <li>Medical</li> <li>Healthcare</li> <li>Dental</li> </ul>	<ul style="list-style-type: none"> <li>General Use</li> <li>General Lab Work</li> <li>Life Sciences</li> <li>Chemical</li> </ul>	<ul style="list-style-type: none"> <li>Medical</li> <li>Healthcare</li> <li>Dental</li> </ul>	<ul style="list-style-type: none"> <li>General Use</li> <li>Life Sciences</li> <li>Semiconductor</li> </ul>
---	--	---	---

## STANDARDS & CERTIFICATIONS

----------	--------------	------------------	--

## COLOUR

<p>Natural White</p> <p>Green (Aloe Vera)</p>	<p>Light Yellow</p>	<p>Light Yellow</p>	<p>Yellow</p>
---	---------------------	---------------------	---------------

## WEIGHT

5.0g   6.0g	5.5g   6.5g   7.0g   8.0g	6.0g   6.3g   6.5g	5.5g   6.0g   7.0g   8.0g
-------------	---------------------------	--------------------	---------------------------

## THICKNESS

\*Depends on glove weight

Wall : Min 0.07mm Palm : Min 0.11mm Finger : Min 0.13mm	Wall : Min 0.07mm Palm : Min 0.10mm Finger : Min 0.12mm	Wall : Min 0.08mm Palm : Min 0.10mm Finger : Min 0.12mm	Wall : Min 0.08mm Palm : Min 0.10mm Finger : Min 0.12mm
---	---	---	---

### MATERIAL

Natural Rubber

### AVAILABLE SIZES

XS: 75 ± 5mm | S: 85 ± 5mm | M: 95 ± 5mm | L: 105 ± 5mm | XL: 115 ± 5mm

### LENGTH

Min 240mm

### PACKAGING

100 Gloves/Dispenser | 10 Dispenser/Case | 1,000 Gloves/Case