

NATURAL RUBBER

# POLYMER LATEX GLOVES

These powder-free, biodegradable natural rubber latex gloves provide enhanced tactile sensitivity, excellent comfort, and a secure grip for various tasks and environments. Available in multiple textures, they are ambidextrous with a beaded cuff and comply with FDA food safety standards for reliable use.



Powder Free



Ambidextrous



Polymer Coated  
(Beaded: 240mm)



Weight Range  
5.0g ~ 6.0g

## Standards & Certifications



FDA CFR 177.2600  
- US FOOD CONTACT



## Colour

Natural White

Green (Aloe Vera)

## Application



Medical



Healthcare



Dental



Scan the QR Code to  
Download Our Latest Catalogue

# LATEX GLOVES

Polymer Latex Gloves	Chlorinated Latex Gloves	Chlorinated Latex EXAM Gloves	Double Chlorinated Gloves
<ul style="list-style-type: none"> <li>Complies to FDA food handling requirements</li> </ul>	<ul style="list-style-type: none"> <li>Complies to FDA food handling requirements</li> <li>Good barrier protection against chemicals, virus and dangerous substances</li> </ul>	<ul style="list-style-type: none"> <li>Good barrier protection against chemicals, virus and dangerous substances</li> <li>Complies to EU MDR requirements</li> </ul>	-

- Powder Free
- Latex gloves designed to provide enhanced tactile sensitivity
- Soft, elastic latex material provide excellent comfort
- Available in finger textured, smooth, and textured surface
- Ambidextrous, beaded cuff



## APPLICATIONS

<ul style="list-style-type: none"> <li>Medical</li> <li>Healthcare</li> <li>Dental</li> </ul>	<ul style="list-style-type: none"> <li>General Use</li> <li>General Lab Work</li> <li>Life Sciences</li> <li>Chemical</li> </ul>	<ul style="list-style-type: none"> <li>Medical</li> <li>Healthcare</li> <li>Dental</li> </ul>	<ul style="list-style-type: none"> <li>General Use</li> <li>Life Sciences</li> <li>Semiconductor</li> </ul>
---	--	---	---

## STANDARDS & CERTIFICATIONS

----------	--------------	------------------	--

## COLOUR

<p>Natural White</p> <p>Green (Aloe Vera)</p>	<p>Light Yellow</p>	<p>Light Yellow</p>	<p>Yellow</p>
---	---------------------	---------------------	---------------

## WEIGHT

5.0g   6.0g	5.5g   6.5g   7.0g   8.0g	6.0g   6.3g   6.5g	5.5g   6.0g   7.0g   8.0g
-------------	---------------------------	--------------------	---------------------------

## THICKNESS

\*Depends on glove weight

Wall : Min 0.07mm Palm : Min 0.11mm Finger : Min 0.13mm	Wall : Min 0.07mm Palm : Min 0.10mm Finger : Min 0.12mm	Wall : Min 0.08mm Palm : Min 0.10mm Finger : Min 0.12mm	Wall : Min 0.08mm Palm : Min 0.10mm Finger : Min 0.12mm
---	---	---	---

### MATERIAL

Natural Rubber

### AVAILABLE SIZES

XS: 75 ± 5mm | S: 85 ± 5mm | M: 95 ± 5mm | L: 105 ± 5mm | XL: 115 ± 5mm

### LENGTH

Min 240mm

### PACKAGING

100 Gloves/Dispenser | 10 Dispenser/Case | 1,000 Gloves/Case